Ceka med Gelbildender Faserverband / Gelling Fiber Dressing

Beschreibung / Description

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Ceka med Gelling Fiber Dressing is a soft, conformable, non-woven pad or ribbon dressing composed of hydrocolloid fibers and enhanced fibers. This conformable and highly absorbent dressing absorbs wound fluid and creates a soft gel, maintains a moist environment which supports the body's healing process and aids in the removal of unnecessary material from the wound (autolytic debridement), without damaging newly formed tissue. This primary dressing should be used with a secondary cover dressing.

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Ceka med Gelling Fiber Dressing is used for the management of moderately to heavily exuding acute or chronic wounds. It can also be pre-moistened for effective use on dry or lightly draining wounds, including: Lower leg ulcers, pressure ulcers (Stage II-IV) and diabetic ulcers; surgical wounds (post-operative, donor sites, wounds left to heal by secondary intent, dermatological); Partial thickness burns; traumatic wounds (abrasions and lacerations); oncology wounds (if moderately or heavily exuding, superficial or deep).

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Patients with allergies to the sodium carboxymethyl cellulose (CMC) should refrain from using this product.

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No adverse effects have been reported to date with correct use.

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1. Before applying the dressing, cleanse the wound area with an appropriate wound cleaner.

2. Apply the dressing onto the wound and cover with a moisture retentive dressing, like foam dressing, gauze, or other appropriate dressing. Ceka med Gelling Fiber Dressing should overlap 1,0 cm or 1/2 inch onto the skin surrounding the wound.

3. When using Ceka med ribbon in deep cavity wounds, only fill deep wounds to 80 %, as the ribbon will expand to fill the wound space on contact with wound fluid.

4. In donor sites, Ceka med Gelling Fiber Dressing may be left in place for up to 14 days.

5. All wounds should be inspected frequently. Remove the Ceka med Gelling Fiber Dressing when clinically indicated (i.e., leakage, excessive bleeding, suspicion of infection).

6. The Ceka med Gelling Fiber Dressing is designed to remain in place not to exceed 7 days. The dressing should be changed when it is saturated with wound fluid or if the cover dressing is leaking or the cover dressing's edges are bunching or rolling up.

7. When using Ceka med Gelling Fiber Dressing to wounds with lower exudate levels, moisten it with sterile physiological saline solution (0.9 %) before changing the dressing. Any gel residue on the wound should be removed when cleansing the wound. Deep wounds in particular should be well irrigated.

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In case of wound swelling or redness for more than 48 hours after use of Ceka med Gelling Fiber Dressing, contact a doctor.

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1. Do not reuse. Do not use the dressing if the sterilized packaging is damaged or open.

2. This product is single use only and is supplied sterile. Discard any unused portion of the product after dressing the wound.

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1. Keep in a clean and dry place in specified temperature range and away from sunlight.

2. Disposal should be handled according to local environmental procedures.

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5 years from date of manufacture