Ceka med Silikonschaumverband / *Silicone Foam Dressing Non-Border*

Beschreibung / Description

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Ceka med Silicone Foam Dressing Non-Border has triple-layered structure which consists of a soft silicone wound contact layer, a absorbent polyurethane foam layer and a vapour permeable and waterproof polyurethane film layer. Ceka med Silicone Foam Dressing Non-Border can absorb exudates, maintain a moist wound environment and minimise the risk for maceration. Dressing material content: Silicone and polyurethane.

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Ceka med Silicone Foam Dressing Non-Border is indicated for a wide range of exuding wounds such as pressure ulcers, leg and foot ulcers, traumatic wounds, and skin tears.

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1. It should not be used on persons with a known allergy or hypersensitivity to the dressing or its components.

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Ceka med Silicone Foam Dressing Non-Border may be used by lay person under supervision of health care professionals.

- 1. Cleanse the wound according to clinical practice. Dry the surrounding skin thoroughly.*
- 2. Select an appropriate size/shape. The wound pad should cover the surrounding skin by at least 1-2 cm.*
- 3. Remove the first release liner and apply the adherent side to the wound.*
- 4. Remove the remaining release liner(s) and smooth down the border on the skin. Do not stretch the dressing.*
- 5. The dressing change interval may be several days, or change the dressing before it is fully saturated, at signs of leakage, or as indicated by clinical practice. When necessary, fixate the dressing with bandage, tapes or compression.*

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- 1. Do not use the product together with oxidizing agents such as hypochlorite solutions or hydrogen peroxide.*
- 2. If reddening or sensitization occur, discontinue use and consult a health care professional for appropriate treatment.*
- 3. Do not reuse. If reused, performance of the product may not meet expectations and cross contamination may occur.*
- 4. Sterile. Do not use if sterile packaging is damaged or opened prior to use. Do not re-sterilize.*
- 5. In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary on these population groups, this dressing should be used with caution following a clinician's recommendation.*

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- 1. Keep in a clean and dry place in specified temperature range and away from sunlight.*
- 2. Disposal should be handled according to local environmental procedures.*

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3 years from date of manufacture