## Ceka med Silikonschaumverband Border / *Silicone Foam Dressing Border*

### Beschreibung / *Description*

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*Ceka med Silicone Foam Dressing Border has an advanced five-layer structure which consists of a soft silicone wound contact layer, a super absorbent polyacrylate fiber layer, a non-woven distributing layer, a polyurethane foam layer, and a vapour permeable and waterproof polyurethane film layer. Ceka med Silicone Foam Dressing Border is highly absorbent dressing that absorbs exudates and maintains a moist wound environment. Dressing material content: Silicone, polyurethane, polyacrylate, fluff pulp and polyolefin.*

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*Ceka med Silicone Foam Dressing Border is designed for a wide range of exuding wounds such as pressure ulcers, leg and foot ulcers, traumatic wounds, surgical wounds and skin tears. The dressing may also be used as part of a prophylactic therapy to help prevent pressure ulcers.*

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*1. It should not be used on persons with a known allergy or hypersensitivity to the dressing or its components.*

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*Ceka med Silicone Foam Dressing Non-Border may be used by lay person under supervision of health care professionals.*

- 1. Cleanse the wound according to clinical practice. Dry the surrounding skin thoroughly.*
- 2. Select an appropriate size/shape. The wound pad should cover the surrounding skin by at least 1-2 cm.*
- 3. Remove the first release liner and apply the adherent side to the wound.*
- 4. Remove the remaining release liner(s) and smooth down the border on the skin. Do not stretch the dressing.*
- 5. The dressing change interval may be several days, or change the dressing before it is fully saturated, at signs of leakage, or as indicated by clinical practice. When necessary, fixate the dressing with bandage, tapes or compression.*

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- 1. Do not use the product together with oxidizing agents such as hypochlorite solutions or hydrogen peroxide.*
- 2. If reddening or sensitization occur, discontinue use and consult a health care professional for appropriate treatment.*
- 3. Do not reuse. If reused, performance of the product may not meet expectations and cross contamination may occur.*
- 4. Sterile. Do not use if sterile packaging is damaged or opened prior to use. Do not re-sterilize.*
- 5. In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary on these population groups, this dressing should be used with caution following a clinician's recommendation.*

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